


DGC-600W


Item code : 100001216005



DGC-600W Clamp Holder
Item code : 100001216019

2 holders necessary to hold
one degassing chamber

2 holders in 1 pack



For Hi-Flow. Large capacity.

【Specification】

Maximum Flow Rate : 600 mL/min

(O₂ removal >50%, 16.6kPa vacuum, water 25°C)

Pressure Drop: <20kPa (at 600mL/min water 25°C)

Fluid type: Deionized water, pure water (neutral buffer solvent)

Wetted Materials: Silicone, PVC, NBR, Polyamide

External Dimensions: 4.49in (W) x 111in (L)

Fluidics Port connection: Inlet: 10mm Barb Outlet: 10mm Barb

Vacuum Port connection: Barb connector for 3mm ID elastomeric tubing

Pressure resistance: 0.1 MPa

Operating temperature: 5-40 (°C)

